

CD-1 Color Detector

Detect and sort specially marked animals

Phason's CD-1 Color Detector detects painted animals and sorts them according to special settings.

The advanced-technology sensor mounts in the sorter and connects quickly and easily to the auxiliary inputs on the PS-1 Scale Head.



Uses for color detection



Uses for the CD-1 include gilt selection, sick pig segregation, and more.

For example, you can select and mark gilts for breeding. When one of the painted gilts moves through the sorter, the CD-1 detects the painted animal, and the system sorts it according to the color sorting settings.

This allows you not only to separate the gilts, but also to confirm that they are the correct weight.

Features

- ◆ Quick and easy installation in new or existing Pro-Sort scales
- ◆ Easy configuration using Phason's OMNI software
- ◆ 90-day limited warranty

Additional information

For more information about Pro-Sort, OMNI, or other Phason products, contact your dealer or visit www.phason.ca.